



CHANTARELLE
Laboratory Derm Aesthetics



THERMOLED MASTER®

24-programme non-invasive device
for integrated therapies



INNOVATIVE TREATMENTS

THERMO method - thermo-modelling

photodynamic method PDT LED:

IR 850nm, R 650nm, Green 525nm, UV Blue 420nm

THERMOLED MASTER

IMMEDIATE EFFECTS + LONG-TERM BIO-STIMULATION

24-programme non-invasive device for integrated therapies:

- using the THERMO method - thermo-modelling and restructuring skin collagen fibres with rotating and pulsating magnetic field; 39 degrees C, 41 degrees C, 42 degrees C

WITH DERMO-THERAPY FOR THE CONCOMITANT SKIN CONDITION

- using the photodynamic method PDT LED: IR 850nm, R 650nm, Green 525nm, UV Blue 420nm

TREATMENTS OFFERED

WITHIN INTEGRATED PROGRAMMES:

- modelling and lifting, anti-ageing, anti-deep wrinkles,
- for dry, rough skin, with symptoms of irritation,
- for oily, problematic skin, with symptoms of inflammation,
- for capillary, sensitive and reddened skin,
- for skin with discolouration, uneven skin tone and scars,
- special programmes for under-eye skin

24 READY-TO-USE PROGRAMMES:

for face skin problems, for skin around eyes and for local body areas:

18 programmes THERMO + PDT LED

+

6 programmes PDT LED



THERMOLED MASTER combines 5 technologies:
THERMAL WITH ROTATING AND PULSATING MAGNETIC FIELD LED LIGHT
COMBINED THERMAL
+ LED LIGHT
+ PHOTODYNAMIC PDT LED
WITH FORMULAS CONTAINING CHANTARELLE BIO-CHROMOPHORES
THERMAL METHOD WITH DEDICATED CHANTARELLE THERMO PRODUCTS

In the device, you can disable the THERMO method - rotating and pulsating magnetic field depending on potential contraindications.

SPECTACULAR EFFECTS OF TREATMENTS USING THERMOLED MASTER:

- face oval modelling/lifting
- smoothed wrinkles
- local body area modelling
- reduction of discolouration
- evened out skin tone
- moisturisation, reduction of roughness and symptoms of irritation
- sebum regulation and reduction of symptoms of inflammation
- reduction of swelling and puffiness



THERMOLED MASTER PROGRAMS

Head/thermal plate

Progr. no.	THERMO programs	Program duration	Magnetic field parameters	Stabilised temperature
P1.	THERMO switched off	20 mins	-----	-----
P2.	THERMO-INFUSION /WOMEN&MEN/ FACE	15 mins	Modulated rotating magnetic field with controlled rotation speed + Pulsating magnetic field	39°C Temperature quickly felt on the plate
P3.	THERMO-LIFTING LEVEL 1 /WOMEN&MEN/ FACE&BODY	30 mins		41°C Temperature quickly felt on the plate
P4.	THERMO-LIFTING LEVEL 2 FORTE /WOMEN&MEN/ BODY	40 mins		42°C Temperature quickly felt on the plate

Head/ LED ring

Progr. no.	Description of the program, indications for use	Duration	IR	R	G	B
			850nm	650nm	525nm	420nm
P.1	Deep furrows and wrinkles, moisturisation and restoration of the epidermal barrier, reduction of the appearance of scars	LED program time equal to the thermal program - automatic adjustment	X	X		
P.2	Anti-wrinkle action, sebum normalisation and hydration, evening out of skin tone, increase in skin density			X	X	X
P.3	Visible capillaries, redness, spots and discolourations, moisturisation, smoothing				X	X
P.4	Keratinisation disorders, roughness and dryness of the epidermis; strong hydration, strengthening of the barrier; under-eye program				X	X
P.5	Oily, problematic, dehydrated skin with active inflammatory conditions		X		X	X
P.6	Comprehensive anti-ageing program Face, neck, décolletage and local body areas.		X	X	X	X

DEVICE ADVANTAGES:

- Innovative technologies:
 - THERMO thermal method
 - magneto-phototherapy LED
 - pulsating magnetic field
 - rotating magnetic field
 - targeted photodynamic therapy PDT
- A very broad treatment scope (24 programmes) – increasing the attractiveness of the offering of the beauty salon
- Immediate effects – visible immediately after the first treatment
- Long-term bio-stimulation effect
- Non-invasive treatment procedures – NO SCALPEL, NO PAIN, NO STRESS
- The patient may return to daily routines immediately after the treatment
- Safety and comfort – for the patient and for the beauty salon, minimisation of side effects compared to invasive treatments
- Quick return on investment
- No need of special maintenance
- Simple and easy-to-follow treatment procedures
- Intuitive interface and pre-programmed treatments
- Option to combine treatments in combined therapies – this broadens the offering of the beauty salon





THERMOACTIVE PREPARATIONS designed for THERMOLED MASTER®:

Head - work in thermo mode:

- ◆ **ABSOLUTE RICH MOISTURE**
PDT THERMO ABSOLUTE MOISTURISING CERAMIDE SERUM Hydrolipidic Firmness Booster (CP184330)
PDT THERMO ABSOLUTE MOISTURISING LIFT EFFECT MASK Ceramide & Collagen Booster (CP0861)
- ◆ **SUPERIOR YOUTH BOOST**
PDT NIACINAMIDE & GABA Boost-Serum Recovery Wrinkle Filler (CP184130)
PDT PERFECT SMOOTH RICH RECOVERY BOOST-MASK with niacinamide (CP0859)
- ◆ **SEBUMELAN HOLISTIC**
THERMO CITRIC ENZYME PEEL pH 4.5, Normalising Anti Acne Spot & Scars (CP0420)
PDT THERMO ANTI-BACTERIAL SERUM pH 4.5, Normalising Anti Acne Spot & Scars (CP1895)
- ◆ **AGEMELAN HOLISTIC**
PDT THERMO FACE & NECK SERUM pH 4.5 (CP1891)
PDT THERMO DECOLETTE & BODY SERUM pH 4.5 (CP1892)
PDT THERMO EYE SERUM pH 4.5 (CP0637)
- ◆ **THERMO-LIFTING**
THERMO-LIFTING KIT Face & Body Skin (CP1511)
- ◆ **THERMO-MODELLING**
SLIM-ESTHETIC CINNAMON CREAM THERMO-ACTIVE (CP1881)
CELLU-ESTHETIC GRAPEFRUIT CREAM THERMO-ACTIVE (CP1883)

Combined treatments:

after applying the PDT laser or a manual photodynamic preparation, wait until it is absorbed, then apply the THERMO-LIFTING kit with the THERMOLED MASTER® device

Applies to photodynamic preparations from the lines:

BLACK ESSENTIAL

LIFT PEPTANGO R (CP1145, CP0540)

