

# THERMOLED MASTER<sup>®</sup>

24-programme non-invasive device  
for integrated therapies



THERMO method - Thermo-modelling  
Photodynamic method  
PDT LED:  
IR 850nm, R 650nm,  
Green 525nm, UV Blue 420nm



# THERMOLED MASTER

IMMEDIATE EFFECTS + LONG-TERM BIO-STIMULATION

## 24-programme non-invasive device for integrated therapies:

- using the THERMO method - thermo-modelling and restructuring skin collagen fibres with rotating and pulsating magnetic field; 39 degrees C, 41 degrees C, 42 degrees C

### WITH DERMO-THERAPY FOR THE CONCOMITANT SKIN CONDITION

- using the photodynamic method PDT LED: IR 850nm, R 650nm, Green 525nm, UV Blue 420nm

### TREATMENTS OFFERED

#### WITHIN INTEGRATED PROGRAMMES:

- modelling and lifting, anti-ageing, anti-deep wrinkles,
- for dry, rough skin, with symptoms of irritation,
- for oily, problematic skin, with symptoms of inflammation,
- for capillary, sensitive and reddened skin,
- for skin with discolouration, uneven skin tone and scars,
- special programmes for under-eye skin

#### 24 READY-TO-USE PROGRAMMES:

for face skin problems, for skin around eyes  
and for local body areas:

18 programmes THERMO + PDT LED

+

6 programmes PDT LED

THERMOLED MASTER combines 5 technologies:  
THERMAL WITH ROTATING AND PULSATING MAGNETIC FIELD LED LIGHT  
COMBINED THERMAL  
+ LED LIGHT  
+ PHOTODYNAMIC PDT LED  
WITH FORMULAS CONTAINING CHANTARELLE BIO-CHROMOPHORES  
THERMAL METHOD WITH DEDICATED CHANTARELLE THERMO PRODUCTS

In the device, you can disable the THERMO method - rotating and pulsating magnetic field depending on potential contraindications.

### SPECTACULAR EFFECTS OF TREATMENTS USING THERMOLED MASTER:

- face oval modelling/lifting
- smoothened wrinkles
- local body area modelling
- reduction of discolouration
- evened out skin tone
- moisturisation, reduction of roughness and symptoms of irritation
- sebum regulation and reduction of symptoms of inflammation
- reduction of swelling and puffiness







## THERMOLED MASTER PROGRAMS

### Head/thermal plate


Progr. no.	THERMO programs	Program duration	Magnetic field parameters	Stabilised temperature
P1.	THERMO switched off	20 mins		
P2.	THERMO-INFUSION /WOMEN&MEN/ FACE	15 mins	Modulated rotating magnetic field with controlled rotation speed + Pulsating magnetic field	39°C Temperature quickly felt on the plate
P3.	THERMO-LIFTING LEVEL 1 /WOMEN&MEN/ FACE&BODY	30 mins		41°C Temperature quickly felt on the plate
P4.	THERMO-LIFTING LEVEL 2 FORTE /WOMEN&MEN/ BODY	40 mins		42°C Temperature quickly felt on the plate

### Head/ LED ring

Progr. no.	Description of the program, indications for use	Duration	IR 850nm	R 650nm	G 525nm	B 420nm
P.1	Deep furrows and wrinkles, moisturisation and restoration of the epidermal barrier, reduction of the appearance of scars	LED program time equal to the thermal program - automatic adjustment	X	X		
P.2	Anti-wrinkle action, sebum normalisation and hydration, evening out of skin tone, increase in skin density			X	X	X
P.3	Visible capillaries, redness, spots and discolourations, moisturisation, smoothing				X	X
P.4	Keratinisation disorders, roughness and dryness of the epidermis; strong hydration, strengthening of the barrier; under-eye program			X	X	
P.5	Oily, problematic, dehydrated skin with active inflammatory conditions		X		X	X
P.6	Comprehensive anti-ageing program Face, neck, décolletage and local body areas		X	X	X	X

## DEVICE ADVANTAGES:

- Innovative technologies:
  - THERMO thermal method
  - magneto-phototherapy LED
  - pulsating magnetic field
  - rotating magnetic field
  - targeted photodynamic therapy PDT
- A very broad treatment scope (24 programmes) – increasing the attractiveness of the offering of the beauty salon
- Immediate effects – visible immediately after the first treatment
- Long-term bio-stimulation effect
- Non-invasive treatment procedures – NO SCALPEL, NO PAIN, NO STRESS
- The patient may return to daily routines immediately after the treatment
- Safety and comfort – for the patient and for the beauty salon, minimisation of side effects compared to invasive treatments
- Quick return on investment
- No need of special maintenance
- Simple and easy-to-follow treatment procedures
- Intuitive interface and pre-programmed treatments
- Option to combine treatments in combined therapies – this broadens the offering of the beauty salon







## THERMOACTIVE PREPARATIONS designed for THERMOLED MASTER®:

### Head - work in thermo mode:

#### ◆ ABSOLUTE RICH MOISTURE

PDT THERMO ABSOLUTE MOISTURISING CERAMIDE SERUM Hydrolipidic Firmness Booster (CP184330)  
PDT THERMO ABSOLUTE MOISTURISING LIFT EFFECT MASK Ceramide & Collagen Booster (CP0861)

#### ◆ SUPERIOR YOUTH BOOST

PDT NIACINAMIDE & GABA Boost-Serum Recovery Wrinkle Filler (CP184130)  
PDT PERFECT SMOOTH RICH RECOVERY BOOST-MASK with niacinamide (CP0859)

#### ◆ SEBUMELAN HOLISTIC

THERMO CITRIC ENZYME PEEL pH 4.5, Normalising Anti Acne Spot & Scars (CP0420)  
PDT THERMO ANTI-BACTERIAL SERUM pH 4.5, Normalising Anti Acne Spot & Scars (CP1895)

#### ◆ AGEMELAN HOLISTIC

PDT THERMO FACE & NECK SERUM pH 4.5 (CP1891)  
PDT THERMO DECOLETTE & BODY SERUM pH 4.5 (CP1892)  
PDT THERMO EYE SERUM pH 4.5 (CP0637)

#### ◆ THERMO-LIFTING

THERMO-LIFTING KIT Face & Body Skin (CP1511)

#### ◆ THERMO-MODELLING

SLIM-ESTHETIC CINNAMON CREAM THERMO-ACTIVE (CP1881)  
CELLU-ESTHETIC GRAPEFRUIT CREAM THERMO-ACTIVE (CP1883)

### Combined treatments:

after applying the PDT laser or a manual photodynamic preparation, wait until it is absorbed,  
then apply the THERMO-LIFTING kit with the THERMOLED MASTER® device

Applies to photodynamic preparations from the lines:

BLACK ESSENTIAL

LIFT PEPTANGO R (CP1145, CP0540)



Treatment procedures are provided in the Catalogue of Treatments and Appliances CHANTARELLE®.